



## EMERGENCY EVACUATION TRAINING MANIKIN

### FOR PROFESSIONALS WHO SAVE LIVES

Ruth Lee is a family business whose core values revolve around the care of our employees, service to our customers and uncompromising product quality. We are proud to serve YOU, the dedicated professionals who make such a valuable contribution to society.

We employ highly skilled people and source the most durable materials to ensure our products meet the highest international standards.

Our manikins are designed to represent the correct anatomic weight distribution of a real person, and we offer a range of different models to suit every conceivable scenario. This ensures that your training is as realistic and as effective as possible; preparing you for all eventualities.

LEARN MORE ABOUT OUR  
EXTENSIVE RANGE OF  
TRAINING MANIKINS,  
BESPOKE BAGS,  
OPERATIONAL  
PRODUCTS AND FIRE  
SIMULATION EQUIPMENT

CALL OUR EXPERTS:

+44 (0) 1490 413 282 SALES@RUTHLEE.CO.UK WWW.RUTHLEE.CO.UK



# EMERGENCY EVACUATION TRAINING MANIKIN

## IDEAL FOR THE FOLLOWING TRAINING SCENARIOS:

If there was a genuine emergency in your care home and the people in your care needed to be evacuated from the building, how confident are you that all Personal Emergency Evacuation Plans (PEEPs) would be carried out successfully? Using manikins for your emergency evacuation training lets you try different scenarios so you can be sure that in a real fire, all staff members would be prepared to move someone safely with minimal risk.

## SPECIAL MATERIALS:

Constructed from a hardwearing P200 material, which is both flame retardant and antibacterial - making it ideal for use in care and medical environments.

## FEATURES & BENEFITS:

- This lighter weight manikin gives a sufficient weight to train with but is easily transported in a vehicle if it needs to be shared among different sites.
- Strong webbing loop at the back allows easy handling and storage.
- Soft body joints to eliminate the risk of pinch welts and bruising; common injuries when

using plastic bodied manikins.

- An immobilisation/extrication collar can be fitted around the neck of our adult manikins to promote proper casualty care.
- Manufactured from an antibacterial material which is also flame retardant.
- Supplied with pyjamas which are easily washable for improved hygiene in hospital and care home environments.
- Comes complete with soft boots which are easily replaced if damaged.

## NOT SUITABLE FOR

Use in temperatures above 100°C or flashover scenarios – See our Fire House manikin.

## OPTIONAL EXTRAS:

- Facemask and wig to add realism.

## USEFUL FOR:

- Hospitals and other medical establishments
- Care Homes.
- Specialist Health & Safety Training Organisations.

## AVAILABLE WEIGHTS & SIZES:

DESCRIPTION	MODEL CODE	HEIGHT (M)	WEIGHT (KG)	NATO STOCK #
25 Kg Adult	RLNEE	1.8	25	N/A