



RUTH LEE™



# LIGATURE TRAINING MANIKIN

## FOR PROFESSIONALS WHO SAVE LIVES

Ruth Lee is a family business whose core values revolve around the care of our employees, service to our customers and uncompromising product quality. We are proud to serve YOU, the dedicated professionals who make such a valuable contribution to society.

We employ highly skilled people and source the most durable materials to ensure our products meet the highest international standards.

Our manikins are designed to represent the correct anatomic weight distribution of a real person, and we offer a range of different models to suit every conceivable scenario. This ensures that your training is as realistic and as effective as possible; preparing you for all eventualities.

LEARN MORE ABOUT OUR  
EXTENSIVE RANGE OF  
TRAINING MANIKINS,  
BESPOKE BAGS,  
OPERATIONAL  
PRODUCTS AND FIRE  
SIMULATION EQUIPMENT

CALL OUR EXPERTS:

+44 (0) 1490 413 282

SALES@RUTHLEE.CO.UK

WWW.RUTHLEE.CO.UK



# LIGATURE TRAINING MANIKIN

## IDEAL FOR THE FOLLOWING TRAINING SCENARIOS:

Institutions within Mental Health or Learning Disabilities, Police, Prisons, Care Homes and other high risk environments who face and need to manage risk of suicide and self-harm on a regular basis, must equip their staff with the right training to safely and sensitively manage hanging scenarios when they occur, including the use of a hook knife to properly release a ligature.

## SPECIAL MATERIALS:

Constructed from the same flame retardant material used in Police ballistic/stab vests, offering 4x the strength and abrasion resistance of 16oz canvas or P.V.C. (conforms to OEKO-TEX® Standard 100, BS EN ISO 12947-2 & 14116: 2008) High density polypropylene webbing (25mm rated to 207 Kg, 50mm rated to 440Kg).

## FEATURES & BENEFITS:

- The manikin has anatomical landmarks including the sternum, rib cage, and substernal notch. The trainer also contains teeth, uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, oesophagus, inflatable lungs and stomach. It allows for placement of a nasogastric tube and accepts Combitube®, E.O.A., E.G.T.A., P.T.L. and i-gel® airways.

- Realistic knee joints promoting correct manual handling techniques including preservation of the C-spine.
- Special D-ring feature enables easy suspension without damaging the manikin.

## NOT SUITED FOR:

- Use in hot fire training - See our Fire House manikin.
- Use in water – See our Water Rescue range.

## EXTENSIVELY USED BY:

- Hospitals and NHS Trusts.
- Police & Prison Service.
- Mental Health & Learning Disability Institutions.
- Equivalent organisations across North & South America, Europe, Australia and Asia.

## AVAILABLE WEIGHTS & SIZES:

DESCRIPTION	MODEL CODE	HEIGHT (M)	WEIGHT (KG)	NATO STOCK #
Adult	RLNLAW/TOR	1.8	25	N/A